

Coremus IM Nail Extractor System Key

MANUFACTURER	NAIL	BOLT No.	DRIVERS	Additional info
ACUMED	POLARIS HUMERAL	P8	2.5/3.5MM HEX	
AESULAP	TARGON FEMORAL	P21		
DOWNS	TARGON TIBIAL	P15		
B BRAUN	HUCKSTEP FEM & TIB	P5	3.5MM HEX	
AO ORIGINAL	UNIVERSAL TAPERED	C2 / C4 / C6		
	V THREAD TAPERED	C1 / C3 / C5		
AUSTOFIX	ALL NAILS	P5		
BIOMET (ZIMMERBIOMET)	ANKLE ARTRODESIS, UNIFLEX TIBIAL, HUMERAL AND RETROGRADE FEMORAL PHOENIX ANKLE	P12	3.5.5.0MM HEX T25 (PHOENIX SET SCREW)	
	KNEE ARTHRODESIS, HOLLAND, FEMORAL, UNIFLEX FEMORAL, UNIFLEX RECON	P17	5.0MM HEX	
	PEDIATRIC NAIL	P7	5.0MM HEX	
	PERI TROCHANTERIC NAIL	P16	3.5/5.0MM HEX 6MM SQUARE	
	C1 / C3 / C5			
CITEFFE	FEMORAL & TIBIAL	P15		
CORIN	HALDER HUMERAL	P5		
CREMASCOLI	HEXAGONE TIBIAL	P14	2.5/3.5/5.0MM HEX	
	FEMORAL & TIBIAL	P22	3.5/5.0MM HEX	
DEPUY ACE (ZIMMERBIOMET)	AIM HUMERAL/TIBIAL	P8	2.5MM/ 3.5MM 5.0MM HEX	
	AIM FEMORAL, SUPRACONDYLAR AND A.R.T FEMORAL	P13		
	ATN TROCHANTERIC	CES - SR - 04		
	CHANDLER TIBIAL	P7	N/A	
	VERSANAIL ANKLE FUSION - TTC	P13	3.5MM HEX	
	VERSANAIL TIBIAL & FEMORAL	P13	3.5MM/ 4.5MM/ 5.0MM HEX	
	AFFIXUS	P19	5.0MM HEX / 3.5MM HEX / LAG BOLT EXTRACTOR	REMOVE SET SCREW
ENCORE	TRUE FIX DISTAL FEMORAL, PROXIMAL HUMERAL AND ANKLE FUSION	P8	3.5MM HEX	
	FEMORAL AND TIBIAL	P1		
FRY	HEXAGONE TIBIAL	P11		
ORTHOFIX	CENTRONAIL	P5		
	VERONAIL	P9		
RICHARDS	TIBIAL	P11		
	RECON	P23		
	FEMORAL	P17		
	IM HIP SCREW	P26		

Coremus IM Nail Extractor System Key

<b>MANUFACTURER</b>	<b>NAIL</b>	<b>BOLT No.</b>	<b>DRIVERS</b>	<b>Additional info</b>
SMITH AND NEPHEW	RUSSEL-TAYLOR HUMERAL	P7	3.5MM HEX	
	RUSSEL-TAYLOR FEMORAL AND RECON	P17	3/16"HEX / 5/32" HEX	
	RUSSEL TAYLOR TIBIAL	P12		
	RUSSEL TAYLOR DELTA TIBIAL	P12		
	RUSSEL TAYLOR DELTA FEMORAL AND RECON	P17	3.5MM HEX / 5/32"HEX	
	TRIGEN (ALL)	P12	3/16" HEX / 6.0MM HEX	
	ENDER NAIL	H1	N/A	
	INTRAPLANT GLIDING NAIL (FOR MAIN FEMORAL NAIL)	P24		
	INTRAPLANT GLIDING NAIL (FOR THE GAMMA ARM)	P9		
	TRI GEN HINDFOOT	P12	3/16" HEX (4.7MM EQUIVALENT)	
	TRIGEN INTERTAN	P12	3/16" HEX	<b>CES-SR-08 Lag bolt extractor</b>
STRYKER	T2 HUMERAL	P4	2.5MM HEX	
(PREVIOUSLY	GAMMA (ORIGINAL)	P25 / P15	4.0MM HEX / 5.0MM HEX	
HOWMEDICA	GAMMA 3	P15	5/32" HEX / 8.0MM HEX / GAMMA SCREW REMOVER	REMOVE SET SCREW
OSTEONICS)	HUMERAL	P5	3.5MM/2.5MM LONG THIN SHAFT	
	GROSS/KEMPF FEMORAL & TIBIAL	P21	3.5MM / 5.0MM HEX	
	IC NAIL FEMORAL, TIBIAL, AND T-2 TIBIAL, FEMORAL ANTEGRADE/RETROGRADE	P14	3.5MM HEX / 6.0MM HEX / 4.0MM HEX	
	SEIDEL HUMERAL	P3	2.5MM HEX (LONG THIN SHAFT)	
	ALTA HUMERAL	P8	T20 TORX	
	ALTA FEMORAL	P12	T25 TORX	
	SOLID (UNREAMED) HUMERAL	P3	2.5MM HEX	
	TI NAIL	P14		

Coremus IM Nail Extractor System Key

MANUFACTURER	NAIL	BOLT No.	DRIVERS	Additional info
SYNTHES	SOLID (UNREAMED) TIBIAL	P5	2.5MM HEX	
(DEPUYSYNTHES)	SOLID (UNREAMED) FEMORAL, RECON AND SUPRACONDYLAR	P10/P15	2.5MM HEX	
	UNIVERSAL TIBIAL AND FEMORAL	C1 / C2 / C3 / C4 / C5 / C6	3.5MM HEX	
	CANNULATED FEMORAL, RECON, AND SUPRACONDYLAR	P10/P15	3.5MM HEX	
	FLEXIBLE HUMERAL	P2	2.5MM HEX	
	CANNULATED TIBIAL	P10/P15	3.5MM HEX	
	EXPERT FEMORAL	P15	T25/T40 TORX	
	EXPERT TIBIAL	P15	T25/T40 TORX	
	DISTAL FEMORAL	P15	3.5MM HEX / SCREW REMOVAL SPIRAL BLADE	
	PFN/A (PROXIMAL FEMORAL NAIL)	P15	SCREW REMOVER (PFNA)	T25/Hex 3.5/ T40
	M/DN/ALL	P13	3.5MM HEX	
	SPIRAL BLADE	SPIRAL BLADE OR P10		P10 ON SLAP ATTACHMENT
ZIMMER	ZMS TIBIAL & HUMERAL	P8	5.0MM HEX	
(NOW ZIMMER BIOMET)	ZMS FEMORAL	P18	5.0MM HEX	
	HANSEN STREET, HARRIS, TIBIA, SCHNEIDER	P8	2.5MM HEX	
	NEFF NAIL	H2	N/A	
	ZICKEL	P12	5/32" HEX / 3/16" HEX	
	Natural	P12 - P13	3.5mm HEX -Screws	5.0MM HEX - END CAP
NUVASIVE	PRECISE	P9	3.5MM HEX	Smooth Screws
UNKNOWN	PUGH HIP SCREW	P5		
	VECTIS (GAMMA TYPE)	P17		
	MARCHETTI FEMORAL	P21		
	FRANKFURT	P25		
	FRANKFURT 12MM	P16		
	ISO KUNTSCHNER SLOT	H1/H2		
	BIOMANTI SRL CLOS NAIL	P10		
	ULTIMAX TIBIAL & FEMORAL	P8		
	FIXANO DSS	P20		
	BROKER-WILLIS FEMORAL		4.5MM HEX	
	BROKER-WILLIS TIBIAL		3.5MM HEX	
	PATHFINDER		4.5MM HEX	
	VECSEI SCREW WITH GROSS/KEMPGF NAILS		4.0MM / 8.0MM HEX	

Coremus IM Nail Extractor System Key

ORTHOSOLUTIONS	OXBRIDGE HINDFOOT	P5	HEX	
Integra (now S+N)	Panta	P2		
DEDICATED PROXIMAL FEMORAL LAG SCREW REMOVERS	Screw Remover - For Synthes PFN Femoral Neck Screw	CES-SR-01		
	Screw Remover - For Synthes UFN & CTN Spiral Blade	CES-SR-02		
	Screw Remover - For Stryker / Howmedica Gamma & Gamma3 Lag Screw	CES-SR-03		
	Screw Remover- For DePuy ATN Lag Screw	CES-SR-04		
	Synthes PFNA Blade Extractor	CES-SR-05		
	Screw Remover for Biomet PTN Lag	CES-SR-06		
	Screw remover for Biomet AFFIXUS lag screw	CES-SR-07		
	Screw remover for S & N INTERTAN lag screw	CES-SR-08		